



Patient Specific  
CRANIOPLASTY IMPLANTS



First time made in India by  
G. Surgiwear Limited  
Rasoolpur Jahanganj, Near Hathoda Chauraha  
Shahjahanpur-242001 UP India

Ahmedabad- 07927545219, Delhi 01165454565  
Mumbai-09821232502, Lucknow-09839626644  
Kolkatta-09830367253, Chennai-09840281838

**SURGIWEAR**®

Patient Specific

# CRANIOPLASTY IMPLANTS



for repair of parietal region



In many patients bony part of the skull is either removed or lost in an accident. their head or face becomes disfigured. It becomes difficult for them to live normal life due to disfigurement. They become depressed and withdrawn from normal life.

To help these patients SURGIWEAR brought patient specific cranioplasty implants. These implants are made exactly to fit the patient & give near perfect contour to body of patient. The implants are based on CT scan data of the patient.

The process starts with CT scan data of the patient.


The operating surgeon has to tell location of screw holes in a simple hand drawing.

SURGIWEAR will prepare a drawing / picture of the proposed implant and an estimate of the cost involved. On acceptance of the drawing and cost the implant will be prepared.

Depending upon the size & complexity the implant may take few hours to few days to build.

The implant can be sent by courier as well as personally to the destination.

Material of Construction :  
TITANIUM Ti6Al4V



for repair of temporo-parietal region

