

G-DRESS

g G DRESS

Material: Non woven spun lace 100% Polyester 48 GSM

Non stick wound pad: 100% PP base non woven pad laminated with polymer mesh 180GSM


PRESENTATION

Code	Size	Pad
GD5	6X5 cm	4x3 cm
GD10	6.5x8.5 cm	4x6 cm
GD15	7x13.5 cm	4x11 cm
GD20	9x20 cm	4.5x11 cm
GD25	9x25 cm	4.5x21 cm
GD30	9x30 cm	4.5x26 cm
GD100	9x100 cm	4.5x100 cm

g G DRESS SWIMPROOF

Material: Non woven PP SMS 25 GSM

Non stick wound pad: 100% PP base non woven pad laminated with polymer mesh 180GSM



PRESENTATION

Code	Size	Pad
GD5SP	6X5 cm	4x3 cm
GD10SP	6.5x8.5 cm	4x6 cm
GD15SP	7x13.5 cm	4x11 cm
GD20SP	9x20 cm	4.5x11 cm
GD25SP	9x25 cm	4.5x21 cm
GD30SP	9x30 cm	4.5x26 cm
GD100SP	9x100 cm	4.5x100 cm

be all free from disease & sufferings

SURGIWEAR®
EN ISO 13485: 2012 COMPANY

SURGIWEAR®
EN ISO 13485: 2012 COMPANY

www.surgiwear.co.in

Made in India by:

G. Surgiwear Limited

Rasoolpur Jahanganj, Shahjahanpur-242001

Tel: +91 5842-223818 Fax: +91 5842-222190

e.mail: surgiwear@hotmail.com

Delhi: 65454565, Mumbai: 66994258, Ahmedabad: 27545218